

# ahn Manual Single Channel Pipettes, Variable Volume.

4 free choice pipettes

1 stand 6 places



\*Made in Germany

### Colour Coding



### Ergonomic Shape



**Sterile handling** - fully autoclavable pipette (by steam at 121°C / 15 min)

**Meets following standards:** DIN EN ISO 8655-1 | DIN EN ISO 8655-2 | DIN EN ISO 8655-6

**Warranty:** 2 years

**CE – mark certified**

Ordering Information valid until : **31 / 12 / 2023**

Volume	Test volume	Increments	Inaccuracy	Imprecision	Cat. No.
0.1-2.5 µL	1.25 µL	0.01 µL	3.00 %	3.00 %	8-100-20-9
	2.5 µL		2.50 %	1.60 %	
0.5-10 µL	5 µL	0.10 µL	2.00%	1.00%	8-101-20-9
	10 µL		1.00%	0.80%	
2-20µL	10 µL	0.10 µL	1.50%	1.00%	8-102-20-9
	20 µL		0.90%	0.40%	
5-50 µL	25 µL	0.50 µL	0.80%	0.40%	8-103-20-9
	50 µL		0.60%	0.30%	
10-100 µL	50 µL	0.50 µL	1.00%	0.50%	8-104-20-9
	100 µL		0.80%	0.15%	
20-200 µL	100 µL	1.00 µL	0.70%	0.30%	8-105-20-9
	200 µL		0.60%	0.15%	
50-200 µL	100 µL	1.00 µL	0.80%	0.20%	8-105-21-9
	200 µL		0.60%	0.15%	
100-1000 µL	500 µL	5.00 µL	1.00%	0.40%	8-106-20-9
	1000 µL		0.60%	0.20%	
200-1000 µL	500 µL	5.00 µL	0.75%	0.25%	8-106-21-9

\*Price w/o VAT

Pipette Starter Kit : 4 free choice pipettes + 6-place Pipette Stand (ref. : **8-SK2-99-9**)

# ahn pipet4u

## Pipette Controller



185.00 €

**Contamination control** – replaceable Membrane filter to minimise risk of contamination.

**Speed setting** – adjustment of pipetting Speed on sample viscosity with 3 selectable speeds

**Warranty** – 2 years

**Sterile handling** – autoclavable pipette cone (120°C / 15 min)

**CE mark certified** – CE IVDR

**UV resistant**



Autoclavable tip cone



Ergonomic push button

Ordering Information valid until : **31 / 12 / 2023**

Description	Sales Unit	Cat. No.
pipet4u® Oasis Pipette Controller 0.1-100 mL	Carton/1 pc.	8-300-00-9

We also supply a large range of high-quality AHN Tips at very competitive prices.

# ahn myTip

## Pipette

**Version** – available in bag, rack, and reload system, with filter or without filter, low retention or non low retention...

**Easy identification of pipette tips** – colour-coded rack.

**High quality material** – tips made of autoclave, non-wettable, virgin polypropylene.

**Metal free**

**CE mark certified** – CE IVDR **(made in Germany)**

**Maintain sample integrity** – Sterile version available (beta irradiation).

**Compatibility list** – available on request

**Catalogue & free samples** – available on request



All prices indicated are without VAT & valid until 31/12/2023

If you have any further questions regarding these products, please do not hesitate to contact us.

[www.aescientific.be](http://www.aescientific.be) +32 2 396 19 10 [info@aescientific.be](mailto:info@aescientific.be)