

Best Prices

Hygromycin B, Solid, High Purity, Cell Culture Tested

Ref. : H-1012 Solid



CAS Number: 31282-04-9
 Molecular Weight: 527.5
 Chemical Formula: $C_{20}H_{37}N_3O_{13}$
 Solubility: Water
 Storage Temperature: 2-8°C
 Solubility: 50 mg/ml H₂O

Price: **156.60 € / 1 G** **570.10 € 5X1 G**

Blasticidin S HCl

Ref. : B-1247



CAS Number: 3513-03-9
 Molecular Weight: 458.9
 Chem. Formula: $C_{17}H_{26}N_8O_5 \cdot HCl$
 Solubility: Water
 Storage Temperature: -20°C
 Purity (HPLC): $\geq 98.0 \%$

Price : **311.60 € / 100 MG**

Puromycin Dihydrochloride

Ref. : P-1033



CAS Number: 58-58-2
 Molecular Weight: 544.4
 Chemical Formula: $C_{22}H_{29}N_7O_5 \cdot 2HCl$
 Solubility: Water
 Storage Temperature: -20°C
 Purity (HPLC): $\geq 98.0 \%$

Price : **66.00 € / 25 MG**
 125.50 € / 100 MG

Dithiothreitol [DTT]

Ref. : C-1029



CAS Number: 3483-12-3
 Molecular Weight: 154.25
 Chemical Formula: $C_4H_{10}O_2S_2$
 Solubility: Water
 Storage Temperature: 2-8°C

Price : **52.00 € / 5 G**
 374.00 € / 100 G

Hygromycin B, Solution (50 mg/ml)

Ref. : TX010-960

Appearance: Clear solution
 CAS No. (solid): 31282-04-9
 Storage and shelf life Store: at +2°C to +8°C.
 Hygromycin B is light sensitive. Store protected from light.
 Shipping conditions: Ambient
 Solvent: Water
 Working concentration Recommended final concentration:
 200 – 500 µg/ml

Prices : **84.50 € / 10 ML** **342.00 € / 5x10 ML**

G418 Sulfate Power

Ref. : TT005-920

Appearance : White to off-white powder
 CAS No. : 108321-42-2
 Storage: at +2°C to +8°C.
 Shipping conditions : Ambient

Prices : **116.00 € / 5 G** **503.50 € / 5x5G**



Help needed?

If you have any further questions regarding our products, please do not hesitate to contact us.

Email: info@aescientific.be

Phone: +32 (2) 396 19 10

All prices are without VAT

Prices valid until 31 / 12 / 2026

Order per email: order@aescientific.com

*Offer valid only in Benelux.