

# CappAero96

MULTICHANNEL PIPETTES

**BUY 2 MULTIPIPETTES**

Free of choice

**& GET**

**10% DISCOUNT**



8 Channels price: **425.00 €**

12 Channels price: **465.00 €**

Cat.	n°. Ch.	Volume range	Inaccuracy %	Imprecision %	Color code
C10-8	8	0.5-10µl	5.00/2.00	3.60/1.00	White
C20-8	8	2-20µl	5.00/2.00	3.60/1.00	Grey
C50-8	8	5-50µl	3.20/1.40	1.60/0.60	Red
C100-8	8	10-100µl	2.50/1.30	1.20/0.50	Orange
C200-8	8	20-200µl	1.80/1.20	0.80/0.40	Yellow
C300-8	8	30-300µl	1.66/1.40	0.60/0.30	Black

Cat.	n°. Ch.	Volume range	Inaccuracy %	Imprecision %	Color code
C10-12	12	0.5-10µl	5.00/2.00	3.60/1.00	White
C20-12	12	2-20µl	5.00/2.00	3.60/1.00	Grey
C50-12	12	5-50µl	3.20/1.40	1.60/0.60	Red
C100-12	12	10-100µl	2.50/1.30	1.20/0.50	Orange
C200-12	12	20-200µl	1.80/1.20	0.80/0.40	Yellow
C300-12	12	30-300µl	1.66/1.40	0.60/0.30	Black

\*See also our Combi Stand (2 modules) for up to 3 mechanical pipettes for: 25.00 €

## Help needed?

If you have any further questions regarding this product, please do not hesitate to contact us.  
Email: [info@aescientific.com](mailto:info@aescientific.com)

Phone: +32 (2) 396 19 10

All prices are without VAT

Offer valid until 31 / 12 / 2021 only in Belgium and Luxembourg.

**Order by fax: 0032 (2) 396 06 73 or email: [order@aescientific.com](mailto:order@aescientific.com)**